

Trojan Battery Company

Wheelchair and Scooter Application Guide

Manufacturer	Model	Number of Batteries	Trojan Supergel Model
AM Scooters	All models	2	SG-30
A-BEC	All models	2	SG-30
Adaptive Driving Systems	Model 12	1	SG-50 or SG-70
	Model 14	2	SG-50
Advanced Technology	All models	2	SG-50 or SG-70
Alpha Unlimited	All models	2	SG-50 or SG-70
American Vermerien	All models	2	SG-70
Amigo Mobility International Inc.	Front Drive	1	SG-30
Braun Corporation	Tri-Wheeler	2	SG-50 or SG-70
Bruno Independent Living Aids Inc.	Regal	2	SG-30
	Cub	"	"
	SuperCub Model 34	"	"
	SuperCub Model 44	"	"
	SuperCub Model 46	2	SG-30 or SG-50
Carter Company	Model 824	2	SG-50 or SG-70
Chauffeur Mobility	Chauffeur Scooters	2	SG-30
	Chauffeur Scooters C Series	2	SG-30 or SG-50
	Viva MWD Powerchair	"	"
	Viva Powerchair	2	SG-50 or SG-70
Damaco Inc.	All models	2	SG-30
Dignified Products	Butler	1	SG-30
	Chauffeur	2	SG-30
E. F. Brewer	Portascoot	1	SG-30
Electric Mobility Corporation	Cycle Chair	2	SG-30
	Deluxe Cycle Chair	"	"
	Rascal	"	"
	Rascal Rover	2	SG-50 or SG-70
	Sparky Joystick	"	"
	2-Seat Cycle Chair	"	"
Electric Wheelchair	Feather Life	2	SG-30

Trojan Battery Company

Wheelchair and Scooter Application Guide

Manufacturer	Model	Number of Batteries	Trojan Supergel Model
Everest & Jennings	Carette	2	SG-30
	Hot Wheels	"	"
	Mobie	"	"
	Tempest	"	"
	Kid Power	"	"
	Model 3N	"	"
	Pacer Alpha United	"	"
	14" Growing Marathon	"	"
	14" Growing 3P/3W	"	"
	Sabre	1	SG-50
	Model 34B	"	"
	Traveler Sprint	2	SG-50
	Quest	"	"
	Models 3H, 3P, 3V, 3W	"	"
	Model 32	"	"
	Explorer	2	SG-70
	Lancer	"	"
	Vortex	"	"
	Model 24 Narrow	1	SG-50 or SG-70
Model 33	"	"	
Model 34 Adult	"	"	
Golden Technologies Inc.	Golden Eagle	2	SG-30
	Golden Companion	"	"
Hoveround	All models	2	SG-30
Independence Chair	Sierra	2	SG-30
Invacare Corporation	Dart/Runabout	2	SG-30
	Jaguar	"	"
	Rabbit	"	"
	Power 9000	"	"
	Tri Rolls	"	"
	Tri Scoot I & II	"	"
	Action Arrow	2	SG-50
	Action Ranger	"	"
	Action XT	2	SG-50 or SG-70
	Arrow XT	"	"
Lakematic	All models	1	SG-70
Leisure Lift	Pacesaver- all models	2	SG-30
Medical Resource	Odyssey 600 & 600HC	2	SG-30

Trojan Battery Company

Wheelchair and Scooter Application Guide

Manufacturer	Model	Number of Batteries	Trojan Supergel Model
Medlift & Mobility Inc.	Adventurer	2	SG-30
	Adventurer Plus	"	"
	Adventurer Plus-4	2	SG-70
	Power Chair	2	SG-30
	Power Chair Plus	2	SG-70
Motovator Electric Wheelchair	All models	1	SG-30 or SG-50
Ortho Kinetics Inc.	Bravo	1 or 2	SG-30
	Pony	"	"
	Colt	2	SG-30
	MVP	"	"
	Triumph	"	"
	Explorer	"	"
	Sierra	"	"
	Encore	"	"
	Lark 3 & 4	"	"
	Outsider	2	SG-50 or SG-70
	460 Sierra XL	"	"
Palmer Industries	All models	2	SG-50 or SG-70
Permobil Inc.	Chairman	2	SG-50 or SG-70
Pride Mobility Products Corporation	Partner	1	SG-30
	Celebrity	2	SG-30
	Legend	"	"
	XL Series	"	"
	Sundancer	"	"
	Shuttle	"	"
	Sidekick	"	"
	Victory	"	"
	Rally	"	"
	Dynamo	"	"
	Wrangler	2	SG-90
	Hurricane	"	"
	Jazzy Power Chair Series		
	1100	2	SG-70
	1103 Mini	2	SG-30
	1105 Foldable	2	SG-50
	1115	"	"
1120 Active-Trac	2	SG-70	
1170 XL	"	"	
1420 & 1470	"	"	

Trojan Battery Company

Wheelchair and Scooter Application Guide

Manufacturer	Model	Number of Batteries	Trojan Supergel Model
Quickie Design	BEC 40 Series	2	SG-30
	Standard Series	"	"
	Targa 14"	"	"
	Targa 16" & 18"	2	SG-50
	Zippie P500	"	"
	P100	2	SG-30 or SG-50
	P190	"	"
Ranger All Season Corporation	1x3	4	SG-30
	2x3	"	"
	Solo	2	SG-30
	SoloIV	"	"
	Safari	"	"
	Supreme	"	"
Redman Power Chair	Chief 107SR	2	SG-70
Shoprider USA Inc.	Streamer 888W	2	SG-30
	Sovereign 888-3 & 888-4	"	"
	Sprinter 889-3 & 889-4	"	"
	Streamer 888WS	"	SG-50
	Sprinter Deluxe 889DX2-4 & DX4-5	"	"
Stand Aid	Power Lift	1	SG-30
	Power Drive	2	SG-30
21st Century Scientific	1BB	2	SG-50
	Bounder	2	SG-90
	Bounder Plus	"	"
	Big Bounder	"	"
	1CP - 1FL	"	"
Wheelcare Inc.	Sportster	2	SG-70